



Product Code . DS-M-10148

Mini Mix Vortex Mixer

Description



Mini Mix Vortex Mixer

Description:

Mini Mixer is used for vibration, blend and stirring of sample tissue, cell, bacteria liquid and chemical reagent in life science and physical and chemical analysis field. Multiple tubes are suitable for 0.2~50ml microtubes and tubes or small container with its diameter is within 108mm.

Features:

1. Mini and compact appearance?powerful function?sucker type under-chassis?super shockproof, ideal for high speed operation.
 2. Various shaking block adapters optional.
 3. Multiple mixing modes, touch-operation, continuous operation and timing set.
 4. Shaking block convenient for replacement, stable and reliable, sturdy and durable eccentric bearing design.
 5. Wide speed range:0-3000rpm/min, continuously variable motor and Brush less DC Motor are optional. Stable low speed, powerful high speed.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/