

Product Code . DS-M-10149

## Vertical Rotating Mixer



### Description

---

#### Vertical Rotating Mixer



##### Description

It is suitable for preventing blood clots, latex diagnosis, immunoprecipitation and similar applications. Vertical Rotating Mixer is suitable for a wide range of test tubes in the laboratory, providing efficient and gentle mixing of experiments to keep biological samples in suspension status, such as blood mixing.

##### Parameters:

- Speed Range: 1rpm ~80rpm
  - Timing Range: 1min ~99h59min
  - Dimension: Width 462 x Depth 150 x Height 180mm
  - Net Weight: 2.1kgs
  - Voltage: AC 90 ~230V, 50/60Hz
  - Power: 48 Watt
-

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)