

Product Code . DS-M-10152

## Magnetic Stirrers 3 Litres



### Description

---

#### Magnetic Stirrers 3 Litres

##### Description:

Magnetic stirrers are often chosen in preference to gear-driven motorized magnetic stirrers because they are quieter, more efficient and have no external moving parts.



- Ultra low profile at just 1.5cm tall
- Speed range of 15 to 1,500 rpm
- Spill and splash proof
- Max volume of 800ml aqueous
- Automatic reverse rotation for better mixing
- Max stirrer bar 40mm
- Chemical resistant surface
- Ideal for titration, limited space, field labs and portable

##### Specifications:

--	--	--	--	--

<b>Number of positions</b>	1	1	1	
<b>Plate size</b>	Ø90mm	Ø120mm	Ø135mm	
<b>Max H2O</b>	0.8L	1.5L	3L	
<b>Plate material</b>	PET	PET	Nylon+GF	
<b>Display</b>	None	None	None	
<b>Max Length Stirrer Bar</b>	40mm	40mm	50mm	
<b>Dimensions (WxDxH)</b>	175x120x15mm	150x145x50mm	150x260x80	

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)