



Product Code . DS-M-10153

Magnetic Stirrers DID215

Description

Magnetic Stirrers DID215



Description:

Magnetic stirrers are often chosen in preference to gear-driven motorized magnetic stirrers because they are quieter, more efficient and have no external moving parts.

- Maintenance free brushless DC motor
- Max. stirring quantity of H₂O at 20L
- Speed range of 0-1500 rpm
- Stainless steel work plate cover with ceramic material provides good chemical resistance performance

Specifications:

| | | | |
|------------------------|--------------------------------------|--------------------|--|
| Plate size | 140 x 140mm | 184x184mm (7 inch) | |
| Max H ₂ O | 5L | 10L | |
| Plate material | Aluminium cover with ceramic coating | Glass Ceramic | |
| Display | LCD | None | |
| Max Length Stirrer Bar | 40mm | 80mm | |
| Dimensions WxDxH | 320 x 180 x 108mm | 215x360x112mm | |

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/