



Product Code . DS-RF-10163

Medical Refrigerator-Double Door

Description

Medical Refrigerator-Double Door

Descriptions:

Medical Refrigerator Medical Refrigerator is specialized refrigerating equipment for cold storage of pharmaceuticals as well as biological products and so forth. It is suitable for hospital, drugstore, pharmaceutical factories, sanitation and antiepidemic station and clinics.

Specifications:

External Size(W*D*H) mm		1220*642*1885	1010*750*2033	1225*795*1980	1800*750*1965
Capacity		650L	758L	1000L	1500L
Temp. Range		2~8℃			
Temp. Accuracy		0.1℃			
Control System		Microprocessor Control & Large Screen LED Display			
Alarm		Audible and visual alarm for: High and low temperature, Power failure alarm, System & Sensor failure, Door ajar			
Refrigeration Type		Forced air refrigeration system			
Refrigerant		R134a, CFC Free			
Defrost		Auto defrost			
Material	Construction	Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)			

	Internal material	Cold-rolled steel coated with anti-bacteria powder			SS304
	External material	Cold-rolled steel coated with anti-bacteria powder			SS304
Shelves		Adjustable shelves made of high quality steel wire			
		10 pcs	16 pcs	12 pcs	12 pcs
Consumption		443W	585W	510W	650W
Power Supply		AC220V±10%, 50/60Hz			
		?, 50/60Hz			?, 50/60Hz
Gross weight		188kg	172kg	245kg	190kg

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/