



Product Code . DS-RF-10164

Medical Refrigerator-Single Door

Description

Medical Refrigerator-Single Door

Description:

It is suitable for hospital, drugstore, pharmaceutical factories, sanitation and antiepidemic station and clinics. Medical Refrigerator Medical Refrigerator is specialized refrigerating equipment for cold storage of pharmaceuticals as well as biological products and so forth.

Specifications:

External Size(W*D*H) mm	540*635*1900	660*655*1980
Capacity	260L	360L
Temp. Range	2??8?	
Temp. Accuracy	0.1?	
Control System	Microprocessor Control, Large Screen LED Display	
Alarm	Audible and visual alarm for: High and low temperature, System & sensor failure, Door ajar, Power failure alarm, Low backup battery alarm.	
Refrigeration Type	Forced air refrigeration system	
Refrigerant	R134a, CFC Free	
Condenser & Evaporator	Bundy tube condenser	

Defrost		Frostless design	
Construction	Structure	Unibody design & mono-assembly foaming (Rigid polyurethane insulation material)	
	Internal material	HIPS engineering plastics, abrasion and acid resistant	
	External material	Cold-rolled steel coated with anti-bacteria powder	
Shelf		Adjustable shelves made of high quality steel wire	
		4pcs	5 pcs
Consumption		200W	350W
Power Supply		AC220V±10%, 50/60Hz	
Gross weight(kg)		87kg	116kg

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/