



Product Code . DS-RF-10168

Medical Refrigerator-Single

Description

Medical Refrigerator-Single

Description:

It is suitable for hospital, drugstore, pharmaceutical factories, sanitation and anti epidemic station and clinics. Medical Refrigerator is specialized refrigerating equipment for cold storage of pharmaceuticals as well as biological products and so forth.

Specifications:

External Size(W*D*H) mm	450*445*575	445*470*700	445*470*900	510*515*900	505*525
Internal Size(W*D*H) mm	382*359*415	400*347*528	400*347*708	419*395*738	385*356
Capacity	50L	70L	95L	130L	105L
Temperature Range	2??8?				
Temperature Accuracy	0.1?				
Control System	Microprocessor Control, Large Screen LED Display				
Alarm	Audible and Visual alarm-Standard: High and low Temperature Alarm, Sensor Failure Audible and Visual Alarm-Optional: Power Failure, Door				
Refrigeration Type	Forced air refrigeration system				
Refrigerant	R134a/R600a, CFC Free				
Condenser & Evaporator	High efficiency air-cooled condenser, made of copper				
Defrost	Auto defrost Removing				
Decondensation	Manual				

Door		Foam/Glass			
Construction	Structure	Unibody design & Mono-assembly foaming (Rigid polyurethane insulation)			
	Internal Material	HIPS blister tank			
	External Material	Cold-rolled steel coated with anti-bacteria powder coating			
Consumption		60W	70W	70W	
Power Supply		AC220V±10%, 50/60Hz			
Gross Weight(Kg)		26kg(F)/	43kg(F)/	53kg(F)/	53kg(F)/
		30kg(G)	48kg(G)	58kg(G)	58kg(G)

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/