



Product Code . DS-RF-10172

Blood Bank Refrigerator (Economic Type)

Description

Blood Bank Refrigerator (Economic Type)

Specifications:

External Size(W*D*H)mm		340*230*280	420*240*300	460*300*340	440*270*310	570*310*370	645*395*420	615*415*450
Internal Size(W*D*H)mm		226*140*210	330*170*225	350*220*250	375*215*250	450*225*300	510*300*355	510*315*355
Capacity		6L	12L	17L	20L	33L	54L	58L
Temp. Range		2??8?						
Storage Time (Ambient 30?)		?48h		?24h				?48h
Material	External	PP						
	Internal	PU						
Max. Loading		6kg	8kg	15kg	18kg	25kg	60kg	60kg
Net Weight		1.9kg	1.9kg	3.2kg	3.5kg	3.5kg	7.8kg	8.5kg
Gross Weight		4.9kg	5.7kg	6.5kg	7.5kg	7.5kg	16.2kg	17.5kg

