



Product Code . DS-RF-10173

-25? Freezer-Horizontal Type

Description

-25? Freezer-Horizontal Type

Description:



The Low Temperature Freezer offers a wide variety of research and storage application, such as low temperature scientific experiments, preservation of plasma, bio material, vaccine, biomedical properties of military products.

It is suitable for research institutes, electronic industries, chemical engineering industries, hospitals, sanitation and anti epidemic stations, university laboratories, military industries.

Specifications:

External Size(L*W*H) mm	696*570*850	865*570*850	1098*617*850	1309*705*874	1843*689*848
Capacity	150L	200L	300L	389L	525L
Temperature Range	-10??-25?				
Control System	Microprocessor Control &LED display				
Alarm	High and low temperature; Condenser fault alarm; Power failure alarm; Temperature control fault alarm;			High or low temperature alarm, temperature sensor failure alarm	

Refrigeration Type	Direct refrigeration					
Refrigerant	R600a/R134a, CFC Free			R600a,CFC Free		
Condenser & Evaporator	Internal-mounted condenser and coil type evaporator					
Construction	Whole foaming, freon cyclopentane / isopentane high density polyurethane foam					
	Internal material: Aluminum					
	External material: Cold-rolled steel coated with anti-bacteria powder					
Basket	1 pc	2 pc	2 pc	1 pc		
Consumption	95W	100W	120W	135W	270W	
Power Supply	220V±10%, 50/60Hz					
Gross Weight (kg)	34	37.5	49	68	99	

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/