

Product Code . DS-DO-10182

Constant-Temperature Drying Oven



Description

Constant-Temperature Drying Oven



Features:

- * Stainless steel inner chamber, round angle structure, adjustable shelves.
- * External material: cold-rolled steel with anti-bacteria powder coating.
- * Door with double-layer tempered glass observation window, adjustable air-tight buckle lock, and silicon rubber seal, which provide good seal ability and efficiently prevent heat loss.
- * PID control with LED display, soft touch operational key, timing function and over-temperature protection.

Specifications:

Capacity	30L	70L	140L	270L
Temp. Range	50~200?			
Temp. Fluctuation	±1?			
Ambient Temp.	5~40?			
Timing Range	1~9999min			
Shelves No.	2 pcs		3 pcs	
Power Consumption	850 W	1000 W	2000 W	3000 W

Power Supply	AC110V/220V±10%,50/60Hz		AC220V±10%,50/60Hz	
Internal Size(W*D*H)mm	300*300*350	400*400*450	500*500*550	600*600*750
External Size(W*D*H)mm	630*590*530	680*630*630	830*650*780	930*740*980
Gross Weight(kg)	55	80	110	150

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/