

Product Code . DS-DO-10189

Mini Autoblot Hybridisation Oven



Description

Mini Autoblot Hybridisation Oven

Description:

- Programmable temperature controller for optimal performance
- Variable rotisserie speed: 0 to 80 RPM
- Choice of 1 free case of 35mm bottles
- 4-Position Rotisserie accepts 8 short or 4 long, 35mm ID bottles
- 6-Position Rotisserie accepts 12 short or 6 long, 35mm ID bottles

Technical Specifications:

- Interior Capacity: 0.5 cubic ft.
- Temperature Range: Ambient +5°C to 99°C
- Variation: $\pm 0.5^{\circ}\text{C}$
- Exterior Dimensions: 14 x 17 x 13" (35 x 42 x 32cm)

-
- Interior Dimensions: 7 x 15 x 8" (18 x 37 x 20cm)
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/