



Product Code . DS-I-10194

Infant Incubator

Description

Infant Incubator

Description:

The infant incubator provides newborns with a temperature and humidity adjustable environment similar to that of the mother's womb, which does not apply to the transfer of the baby.

Specifications:

Power Supply	AC220-230V/50Hz,1300VA
Maximum heater output power	350W/240V
Auxiliary power output	AC220-230V/50Hz,3A
Heater instructions	0?100%
Temperature control mode	Box temperature mode; infant mode
Box temperature control temperature range	25?37??37?39?(Leapfrogging Mode)
Skin temperature control temperature range	34?37?37?38??Leapfrogging Mode?
Sensor temperature display range	5?65?

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/