



Product Code . DS-I-10195

Lighting Incubator

Description

Lighting Incubator

Features:

1. Magnetic door sealing.
2. Microprocessor PID control, accurate and stable.
3. CFC-free refrigerant, eco-friendly and highly efficient.
4. The three sides special glasses ensure good heat insulating ability and desired light intensity.

Safety Function:

1. Over temperature protection.
2. Auto-off protection and alarm for abnormal situation.
3. Parameters can be saved automatically when power off or system error.
4. Temperature sensor for open circuit or short circuit protection.

Specifications:

Capacity	150L	250L	300L	400L
Lighting type	2-side illumination			
Temp. range	With illumination: 10~50 °C, Without illumination: 5~50 °C,			
Temp. fluctuation	±1.5°C			
Temp. uniformity	±2°C			
Lighting illumination	0~6000 Lux (6 degree adjustable)			
Ambient temp.	5~35°C			
Continuous working time	≥180h			
Refrigerant	R134a (CFC free, eco-friendly)			
Shelves	3 pcs/adjustable			
Consumption	1000W	1500W	2000W	2500W
Power supply	220V ±10% 50Hz			
Internal size (mm)	420*420*930	470*460*1180	490*475*1280	560*600*1200
External size (mm)	730*580*1570	665*765*1820	680*780*1920	755*905*1840
Gross weight(kg)	148	158	178	200

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/