



Product Code . DS-I-10196

Lighting Incubator (II)

Description

Lighting Incubator (II)

Safety Function:

1. Auto-off protection and alarm for abnormal situation.
2. Temperature sensor for open circuit or short circuit protection.
3. Over temperature protection.
4. Parameters can be saved automatically when power off or system error.

Features:

1. Microprocessor PID control, accurate and stable.
2. CFC-free refrigerant, eco-friendly and highly efficient.
3. Magnetic door sealing.
4. The three sides special glasses ensure good heat insulating ability and desired light intensity.

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  sales@didacscientific.co.uk

 www.didacscientific.co.uk/