



Product Code . DS-I-10197

Climate Incubator

Description

Climate Incubator

Features:

1. Magnetic door sealing.
2. CFC-free refrigerant, eco-friendly and highly efficient.
3. Auto defrost without influence on the inner temperature.
4. The three sides special glasses ensure good heat insulating ability and desired light intensity.
5. Microprocessor PID control, accurate and stable.

Safety Function:

1. Auto-off protection and alarm for abnormal situation.
2. Temperature sensor for open circuit or short circuit protection.
3. Over temperature protection.
4. Parameters can be saved automatically when power off or system error.

Specifications:

Capacity	250L	300L	400L
Lighting Type	3-side illumination		
Temp. Range	With illumination: 10~50 °C, Without illumination: 5~50 °C		
Temp. Fluctuation	±0.3°C		
Temp. Uniformity	±2°C		
Humidity Range	With illumination: 50%~80%, without illumination: 50%~90%		
Humidity Fluctuation	± 3%RH		
Lighting Illumination	0~12000 Lux		
Illumination Control	5 degree adjustable		
Ambient Temp.	5~35°C		
Continuous Working Time	?180h		
Consumption	2500W	3000W	3500W
Power Supply	AC220V ±10%, 50Hz		
Shelves	3 pcs/adjustable		
Internal Size(W*D*H)mm	470*460*1180	490*475*1280	560*600*1200
External Size(W*D*H)mm	880*670*1820	900*660*1920	1000*730*1870
Gross Weight(kg)	148	158	178

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/