



**Product Code . DS-I-10198**

## **Climate Incubator Class-A**

### **Description**

---

#### **Climate Incubator Class-A**

##### **Features:**

1. Microprocessor PID control, accurate and stable.
2. CFC-free refrigerant, eco-friendly and highly efficient.
3. Magnetic door sealing.
4. The three sides special glasses ensure good heat insulating ability and desired light intensity.

##### **Safety Function:**

1. Over temperature protection.
2. Auto-off protection and alarm for abnormal situation.
3. Parameters can be saved automatically when power off or system error.
4. Temperature sensor for open circuit or short circuit protection.

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)