



Product Code . DS-I-10199

Mould Incubator

Description

Mould Incubator



Features:

- * Power on/off switch, more convenient for operation.
- * Double-layer glass observation window, filled with inert gas, to ensure the cooling performance.
- * UV lamp for disinfection, the interlock function for UV lamp and door.
- * High efficiency refrigeration system, fast cooling speed and uniformity temperature inside.
- * Microprocessor temperature control, high accuracy sensor, LED display.
- * Audio and visual alarm for high/low temperature, sensor error, power failure, door ajar, low battery.

Specifications:

Capacity	150L	200L	250L	300L	400L
Temp. Range	0~65°C				
Temp.	±1 °C				
Fluctuation					
Shelves	3 pcs/adjustable				

Consumption	300W	350W	430W	480W	750W
Power Supply	AC110/220V±10%, 50/60Hz				
Standard Accessory	micro printer (for series III), humidifier				
Internal Size(W*D*H)mm	450*420*850	450*450*1000	480*490*1070	520*500*1170	580*540*1270
External Size(W*D*H)mm	590*620*1400	590*640*1550	620*680*1620	660*690*1720	720*740*1820
Gross Weight(kg)	120	138	147	157	170

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/