



Product Code . DS-I-10204

Platelet Incubator-II

Description

Platelet Incubator-II



Features:

- * Microprocessor control system, LCD display. interface for software connection.
- * UV lamp for disinfection, the interlock function for UV lamp and door.
- * Stainless steel inner chamber, round angle structure, easy cleaning.
- * Double-layer glass observation window, filled with inert gas, to ensure the cooling performance.
- * Audio and visual alarm for high/low temperature, stalling speed.
- * High efficiency refrigeration system, fast cooling speed and uniformity temperature inside.

Specifications:

Temperature Control Range	22±2?	
Temperature Fluctuation	±1?	
Oscillation	60±5cycles/minute	
Amplitude	30mm	
Trays	5	7
450ml Blood Bags	Max.10pcs	Max.15pcs

Temp. Control	Microprocessor controller	
Display	LCD	
Number Of Agitator	Single	
Refrigerant	R134a CFC free	
Power Supply	AC220V±10%, 50/60Hz	
Consumption	180W	250W
Internal Size (W*D*H)mm	400*370*550	450*420*850
External Size (W*D*H)mm	540*570*1070	590*620*1370
Net /Gross Weight(kg)	50/85	90/120

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/