



Product Code . DS-I-10205

Platelet Incubators

Description

Platelet Incubator



Features:

- * Power on/off switch, more convenient for operation.
- * Microprocessor temperature control, high accuracy sensor, LED display.
- * UV lamp for disinfection, the interlock function for UV lamp and door.
- * Double-layer glass observation window, filled with inert gas, to ensure the cooling performance.
- * High efficiency refrigeration system, fast cooling speed and uniformity temperature inside.
- * Audio and visual alarm for high/low temperature, sensor error, power failure, door ajar, low battery.

Specifications:

Temperature	22±2?	
Accuracy	±0.1 °C	
Oscillation	60±5cycles/minute	
Amplitude	50mm	
Trays	5	10
450ml Blood Bags	Max.10pcs	Max.20pcs

Temp. Control	Microprocessor controller	
Display	LCD	
Number Of Agitator	Single	
Refrigerant	R134a CFC free	
Power Supply	AC220V±10%, 50/60Hz	
Consumption	101W	180W
Internal Size (W*D*H)mm	427*470*467	607*535*657
External Size (W*D*H)mm	540*605*960	720*670*1150
Net /Gross Weight	98/130kg	110/150kg

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/