



Product Code . DS-WP-10210

Water Purifier (Ultra-Pure Grade)

Description

Water Purifier (Ultra-Pure Grade)

Features:

1. Automatic microcomputer control system and alarm.
2. Large digital display, show all the running state and parameters.
3. Detect the water quality timely.

Specifications:

Water type	Ultrapure water			
Water output	30 L/H	80 L/H	120 L/H	200 L/H
Ultra-pure water quality	Output water quality $\leq 15 \text{ M}\Omega\cdot\text{cm}$			
	TOC (Supply is RO water) 1~10 ppb			
	Bacteria $\leq 1 \text{ CFU/ml}$			
	PH7			
Indications	Power indication,	Power indication,		
	Process indication	Process indication,		

		Machine condition indication		
Automatic Safety	Power-off protection			
And alarm				
	Automatic water-level control			
	Low water pressure protection			
	Auto restart	Visual/audio	Visual/audio	Visual/audio
		alarm	alarm	alarm
Power supply	AC110/220V, 50/60Hz			
External size	340*350*640	405*395*1000	470*470*1130	600*510*1330
(W*L*H)mm				
Gross Weight (Kg)	40	55	85	120

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/