



Product Code . DS-WP-10211

Water Purifier Medium Type (Automatic RO Water)

Description

Water Purifier Medium Type (Automatic RO Water)

Features:

- * RO membrane clean automatically. High/Low pressure protection.
- * Large LED display, show all the running state and parameters.
- * Alarm Function: Water quality, no water, full water, consumables life ending. Automatic diagnosis.

Specifications:

Water Output Type		RO water			
Water Output Speed		45 L/H	63 L/H	94 L/H	125 L/H
Purifying Procedure		PF+AC+AC+RO			
Water Supply Requirement		Tap water: TDS<200 ppm, 5~45°C, 1.0~4.0Kgf/cm ² (if TDS>200ppm, extra pretreatment filter is recommended)			
Pure Water Quality	Ion rejection rate	96%~99% (New RO membrane)			
	Organic rejection rate	>99%, when MW>200 Dalton			
	Particles and bacteria	>99%			

	rejection rate		
	Bacteria	<0.1cfu/ml (with optional 0.2µm PES terminal filter)	
	Particles(>0.2µm)	<1/ml (with optional 0.2µm PES terminal filter)	
Consumption		120W	240W
Power Supply		AC110/220V±10%, 50/60Hz	
External Size (W*D*H)		540*640*1110mm	
Gross Weight		70kg	

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/