



Product Code . DS-BB-10213

Blood Bank Refrigerator

Description

Blood Bank Refrigerator

Description:

- Blood Bank Refrigerators are available in six different sizes and designed to meet strict requirements for storage of whole blood.
- The complete selection of 4°C refrigerators includes the Micro controlled Temperature control having Audio Visual Alarm in case of Temperature Deviation.
- Micro Controlled Temperature controller with Audio Visual Alarm for high & Low Temperature.
- Heavy Duty Stainless Steel (304) Baskets.
- The unit is mounted on castor wheel for easy mobility.
- Extra-strength non CFC refrigeration System
- CFC Free PUF(Poly Urethane Foam) Insulation to ensure temperature stability and save power consumption.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/