



Product Code . DS-BB-10217

Cryo Bath

Description

Cryo Bath

Description:

Cryoprecipitate Bath is specially designed to thaw Frozen Plasma at +4.0°C/ +37.0°C. Inner chamber made of Stainless Steel (304), and outer body of MS duly powder coated with built in heavy duty industrial pump to circulate the water throughout the chamber to get maximum accuracy.

The Micro Controlled temperature controller controls the temperature with accuracy of $\pm 0.2^{\circ}\text{C}$ specially designed Stainless Steel Basket to accommodate Cryo/Plasma Bags.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/