

Product Code . DS-BB-10221

Platelet Incubator and Agitator



Description

Platelet Incubator and Agitator



Description:

Platelet Incubator and Agitator Designed to maintain a stable and uniform temperature of 22.0°C, The user-friendly design includes an acrylic inner door, and see through outer door, simple calibration is possible via the control panel. This, coupled with a visual and audible alarm system, ensures that your samples are both safe and secure, but also easy to access when you need them.

Features:

- Micro Controlled Temperature controller with Audio Visual Alarm for high & Low Temperature and Door Open
- Display of Set and Process value.
- Heavy Duty Stainless Steel Basket.
- The unit is mounted on castor wheel for easy mobility.
- Extra-strength non CFC refrigeration System.
- CFC Free PUF (Poly Urethane Foam) Insulation to ensure temperature stability and save power

consumption.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/