



Product Code . DS-BB-10222

Cryoprecipitate Freezer

Description

Cryoprecipitate Freezer



Description:

Cryoprecipitate Freezer Provides an ideal freezing environment for the preservation of Cryoprecipitate, Plasma, and test samples. The user-friendly design includes a built-in Stainless steel insulated inner door, and simple calibration is possible via the control panel. This unit Designed to maintain a stable and uniform temperature of up to -80.0°C , Vertical and Horizontal with four different sizes. This coupled with a visual and audible alarm system, ensures that your samples are both safe and secure, but also easy to access when you need them.

Features:

- Micro Controlled Temperature controller with Audio Visual Alarm for high & Low
- Display of Set and Process value.
- Heavy Duty Stainless Steel (304) Basket.
- Temperature and Door Open
- The unit is mounted on castor wheel for easy mobility.

-
- Extra-strength non CFC refrigeration System.
 - CFC Free PUF Insulation to ensure temperature stability and save power consumption.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/