

Product Code . DS-BB-10224

Blood Bag Tube Sealer



Description



Every blood transfusion relies on meticulous safety measures. Blood bag tube sealers play a vital role in this process. These small but essential devices create secure seals on the tubes attached to blood bags, preventing leakage and contamination. Fortunately, the UK boasts a range of companies specializing in blood bag tube sealers. Blood bag tube sealer manufacturer UK develop innovative and reliable sealing technology. They prioritize features like fast sealing times and easy operation, ensuring efficient blood processing in hospitals and blood banks. Blood bag tube sealer supplier UK bridge the gap between manufacturers and healthcare facilities. They offer a wide selection of sealers to suit different needs and budgets. Many suppliers also provide excellent after-sales support, ensuring your equipment functions optimally for years to come. By partnering with a trusted blood bag tube sealer exporter UK, you can gain access to this crucial technology on a global scale. Their expertise ensures proper storage and transportation of blood bags, ultimately safeguarding the health of patients worldwide.

Blood Bag Tube Sealer

Description:

- Blood Bag tube Sealer designed to seal the blood bag tube by high frequency dielectric heating system without causing haemolysis, leakage of blood.
- Sealing of each tube takes approx. 1.5 seconds and is thus very efficient.
- The sealing can be achieved in less than 1.5 seconds per tube by simple operation.
- Pulling the sealing section to both sides can readily cut off sealed tube.

-
- Complete sealing of the blood bag tube is annexure by high frequency dielectric heating.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/