

Product Code . DS-BB-10225

Blood Collection Monitor



Description



Blood Collection Monitor

Description:

LCD Blood Collection Monitor is a compact instrument to provide smooth and gentle rocking for homogeneous mixing with anticoagulant without clot formation of blood from a donor.

Special Features:

- Microprocessor based controller system with accuracy of ± 1 ml.
 - Pre-selection of volume / 16x2 Jumbo LCD display with White Back light.
 - Mixing speed of 12 to 14 rpm for gentle and smooth mixing with anticoagulant prevents Clot formation.
 - Meant for stationary and mobile use for easy & safe blood Collection.
 - Audio and visual alarm for High flow (> 180 ml/min) and Low flow (< 20 ml/min) and commencement of collection.
 - End to end gentle blood mixing with coagulant and control of collection time.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/