

**Product Code . DS-BB-10226**

## Single Pan Balance



### Description

---

#### Single Pan Balance

##### Description:

Micro Controlled Blood bank Scale is designed for weighing Blood and Blood Components with display of Weight and Volume.

##### Special Features:



- \* ABC Mounded body.
- \* Auto Calibration.
- \* LED or LCD displays of Weight and Volume with accuracy of  $\pm 1\text{ml/grm}$ .
- \* Micro Controlled based technology for high accuracy.
- \* Tare provision to account for the weight of the blood bag.
- \* Battery support system up to 8 Hrs.

---

**Specification:**

- Capacity : 10 Kg
- Accuracy :  $\pm$  1gm
- Consumption : 15 V A
- Weight : 2.5 Kg approx.
- Input Volt : 230V AC Power
- 6 V 1.2 Ah Lead Sealed Battery

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)