



**Product Code . DS-GB-10227**

## **Glass Beaker**

### **Description**

---

#### **Glass Beaker**

#### **Description:**

Beakers are a great value! They are made of high-quality borosilicate glass that can be heated directly in an open flame. This glassware will withstand typical laboratory thermal variations common to chemistry processes like heating and cooling.

#### **Specifications:**

- o Material: Borosilicate Glass
- o Capacity: 50ml, 100 ml, 150 ml, 250 ml, 400 ml, 500 ml, 600 ml, 1000 ml, 2000 ml
- o Heat Resistant
- o Chemical Resistant
- o Graduated

