



Product Code . DS-GB-10237

Pyrex Beaker, Low Form

Description

Pyrex Beaker, Low Form



Description:

This 30 ml Pyrex low-form Griffin beaker is used in professional labs. It is made from low expansion, Pyrex glass that can be heated directly in an open flame. This beaker features a single graduated metric scale with white enamel graduations and extra-large marking spot. Pyrex beakers for general laboratory use such as mixing, heating, and cooling. The graduations are approximate and not intended for accurate liquid measurement.

Specifications:

- o Material: Pyrex Glass
- o Capacity: 30 ml, 50 ml, 100 ml, 150 ml, 250 ml, 400 ml, 600 ml, 1000 ml, 2000 ml, 3000 ml, 4000 ml
- o Heat Resistant
- o Chemical Resistant
- o Graduated

o Extra-large marking spot

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/