



Product Code . DS-GB-10249

Beaker, Polypropylene

Description

Beaker, Polypropylene, 50 ml

Description:

Polypropylene beakers can handle hot liquids but not an open flame. This Polypropylene beaker is much less expensive than glass and is suitable for many laboratory uses. This beaker is marked with 10 ml graduations and is approximately 6 cm tall and 5 cm in diameter.

Specifications:

- o Material: Polypropylene
- o Chemical Resistant
- o Graduated
- o Capacity: 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml, 2000 ml
- o With Spout

