



Product Code . DS-GB-10255

Plastic Beaker Set of 5

Description



Plastic Beaker Set of 5

Description:

This Set of Plastic Beakers includes five different sizes, ranging from 50 ml to 1000 ml. Plastic Beakers can hold hot liquids, but are not to be used over an open flame or electric hot plate. Plastic beakers are used for mixing. The graduations on these beakers are approximate and not intended for accurate liquid measurement.

Specifications:

- o Material: Polypropylene
- o Capacity: 50ml to 1000 ml
- o Chemical Resistant
- o Graduated
- o With Spout

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/