



Product Code . DS-GB-10257

Pitcher, Polypropylene, 2000 ml

Description



Pitcher, Polypropylene, 2000 ml

Description:

This 2000 ml Polypropylene lab pitcher is similar to a beaker in shape but with a sturdy handle and larger spout for easy pouring. Polypropylene pitchers can handle hot liquids but not an open flame. This pitcher is marked with 50 ml graduations from 50-2000 ml. It is approximately 17 cm tall and 14 cm in diameter.

Specifications:

- o Material: Polypropylene
- o Capacity: 2000 ml
- o Chemical Resistant
- o Graduated
- o With Handle

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/