

**Product Code . DS-GB-10258**

## **Crucible & Beaker Tongs**



### **Description**

---

#### **Crucible & Beaker Tongs**

##### **Description:**

These tongs are useful when handling a hot crucible, evaporating dishes, flasks, and small beakers.

Material: Nickel plated steel.

Length: 8" long.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)