

Product Code . DS-GB-10259

Beaker Tongs



Description

Beaker Tongs

Description:

- o Use Beaker Tongs to handle hot beakers.
- o These nickel-plated steel tongs are 25 cm (10") long and have plastic coated jaws that open up to 11 cm (4.5").
- o They can hold beakers ranging from 100 ml up to 1000 ml.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/