



Product Code . DS-GF-10261

Glass Funnel

Description



Glass Funnel

Description:

This small glass funnel is made from borosilicate glass and features a short stem. Perfect for use in chemistry experiments.

Specifications:

- o Material: Borosilicate Glass
- o Heat Resistant
- o Chemical Resistant
- o Diameter: 75 mm, 90 mm, 120 mm
- o Shape: Conical
- o Used in Laboratory.

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/