



**Product Code . DS-GF-10264**

## Plastic Funnel

### Description

---

#### Plastic Funnel

#### Description:

Two-ounce capacity with 30 mm stem, 3 mm-diameter stem opening. 11 cm filter paper is recommended for this funnel and, when folded, will fit well inside the funnel without overhang.

#### Specifications:

- o Material: Plastic
- o Chemical Resistant
- o Diameter: 35 mm, 65 mm, 100 mm
- o Shape: Conical

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)