



Product Code . DS-GF-10266

Buchner Funnel, 90 mm, Porcelain

Description

Buchner Funnel, 90 mm, Porcelain

Description:

- o A Buchner Funnel is used for vacuum filtration.
- o This porcelain Buchner funnel is 90mm diameter and has a capacity of about 200 ml.
- o The stem tapers from about 22 mm OD at the end to 30 mm OD at the cone.
- o Porcelain is durable, acid and alkali resistant, easy to clean and autoclavable.
- o Be sure to order with 11 cm filter paper, which will need to be cut down.
- o The one-piece design of the porcelain 90 mm Buchner funnel is perfect for vacuum filtrations.
- o The filter platform has dozens of small holes (about 1/8") to provide filtration efficiency.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/