



Product Code . DS-GF-10272

Powder Funnel

Description

Powder Funnel

Description:

These Polypropylene Funnels are best for transferring powdered reagents.

These funnels are steam autoclavable.

The unique feature of these funnels is its parallel stem that minimizes bridging of powder.

Use with 11 cm filter paper.

Specifications:

- o Material: Plastic
- o Diameter: 80 mm
- o Stem opening: 12.5 mm
- o Stem length: 30 mm

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  sales@didacscientific.co.uk

 www.didacscientific.co.uk/