



**Product Code . DS-MC-10274**

## **Glass Measuring Cylinder**

### **Description**

---

#### **Glass Measuring Cylinder**

##### **Description:**

This graduated cylinders are manufactured from borosilicate glass, which withstands typical laboratory temperature variations. The liquid level forms a true meniscus, making them easy to read. They are also chemically-resistant and thermally-resistant for all applications.

##### **Specifications:**

- o Material: Borosilicate Glass
- o Capacity: 10 ml, 25 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml
- o Heat Resistant
- o Chemical Resistant
- o Graduated
- o Used in Laboratory.

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)