



Product Code . DS-GF-10308

Baffled Culture Flask

Description

Baffled Culture Flask

Description:-

- o Manufactured from Borosilicate glass.
- o Ideal for Use of in preparation of culture.
- o Heat Resistance
- o Chemical Resistant
- o Suitable for repeated autoclave cycles at 121°C.
- o The four equidistant baffles provide consistent agitation when used with compatible orbital shakers.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/