



Product Code . DS-GF-10314

Evaporating Flask

Description

Evaporating Flask

Specifications:

- o Material: Borosilicate Glass
 - o Heat Resistant
 - o Chemical Resistant
 - o Graduated
 - o Used in Chemical Laboratory.
 - o Capacity:
 - 100 ml
 - 250 ml
 - 500 ml
 - 1000 ml
 - 5000 ml
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/