



Product Code . DS-GF-10315

Distillation Flask

Description

Distillation Flask

Description:-

o This flask is commonly used in chemical labs to separate liquid mixtures by the process of distillation.

o Material: Borosilicate Glass

o Heat Resistant

o Chemical Resistant

o Capacity: 100mL, Dia: 65mm, Height: 215mm, Side Dia: 7mm, Arm Length:100

o Capacity: 125mL, Dia: 68mm, Height: 215mm, Side Dia: 7mm, Arm Length:100

o Capacity: 250mL, Dia: 85mm, Height: 226mm, Side Dia: 9mm, Arm Length:130

o Capacity: 500mL, Dia: 105mm, Height: 240mm, Side Dia: 9mm, Arm Length:130