



**Product Code . DS-GF-10316**

## **Pear Shaped Flask**

### **Description**

---

#### **Pear Shaped Flask**

##### **Description:-**

It is made of Borosilicate glass for chemical resistance and thermal performance. The pear shaped design allows distillation down to a small bulk which will reduce losses from thermal decomposition. Ideal for semi-micro distillations.

##### **Specifications:**

- o Material: Borosilicate Glass
- o Capacity: 50 ml, 100 ml, 250 ml, 500 ml
- o Heat Resistant
- o Chemical Resistant
- o Graduated
- o Used in Chemical Laboratory

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)