



Product Code . DS-GF-10326

Iodine Flask

Description

Iodine Flask

Description:-

- o These Glass Flask iodine are made from borosilicate glass.
- o All glass Flask iodine come with thick uniform wall thickness and are printed with white enamel.
- o Determination of Iodine with ground joint and stopper.
- o Heat Resistant
- o Chemical Resistant
- o Capacity: 250, Joint Size: 24/29
- o Capacity: 250, Joint Size: 29/32
- o Capacity: 500, Joint Size: 24/29
- o Capacity: 500, Joint Size: 29/32

