



**Product Code . DS-GF-10327**

## **Pear Shaped Flask with Two Neck**

### **Description**

---

#### **Pear Shaped Flask with Two Neck**

#### **Specifications:**

- o Made of Borosilicate Glass.
- o Side neck at an angle of 25°.
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory
- o Capacity: 50mL, Centre Size: 14/23, Side Socket: , Overall Height: 100mm
- o Capacity: 100mL, Centre Size: 14/23, Side Socket: , Overall Height: 100mm

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)