

Product Code . DS-GF-10355

Two Neck Reaction Flask Lid





Description

Two Neck Reaction Flask Lid

Specifications:

- o Material: Borosilicate Glass
- o Two Neck
- o One center and one angled side neck
- o With standard taper 24/40 joints.
- o Heavy wall
- o Fits 150mm flanged vessels.
- o Flange has a flat ground surface allowing it to be greased.
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/