



Product Code . DS-GF-10358

Four Neck Reaction Flask Lid

Description

Four Neck Reaction Flask Lid

Specifications:

- o Lid has a flat ground surface allowing it to be greased or used with an o-ring.
- o Material: Borosilicate Glass
- o Four Neck
- o With standard taper joints.
- o Heavy wall
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/