



**Product Code . DS-GF-10359**

## **Five Neck Reaction Flask Lid**



### **Description**

---

#### **Five Neck Reaction Flask Lid**

#### **Specifications:**

- o Lid has a flat ground surface allowing it to be greased or used with an o-ring.
- o Material: Borosilicate Glass
- o Five Neck
- o Heavy wall
- o Flat Flange and Interchangeable Joint.
- o Low Expansion Borosilicate Glass.
- o Central: 19/26, Parallel: 19/26, Side 5°: 24/29, Side 10°: 19/26, Side 15°: NA
- o Central: 19/26, Parallel: 14/23, Side 5°: 14/23, Side 10°: 14/23, Side 15°: 29/32
- o Central: 29/32, Parallel: 14/23, Side 5°: 14/23, Side 10°: 14/23, Side 15°: 29/32
- o Heat Resistant
- o Chemical Resistant

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)