



Product Code . DS-GF-10362

Flasks-E, Without Side Cut for Dissolution Apparatus

Description

Flasks-E, Without Side Cut for Dissolution Apparatus

Specification:-

- o Smooth hemispherical base with no flaws or defects, straight side and a flat flange
- o Graduated in intervals of 50ml
- o Ideal for determinant the dissolution characterism of medicine in capsules and tablets.
- o Material: Borosilicate Glass
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/